

Here are notes about how to prepare your samples for cell sorting:

1. Cell concentration (before sort) should be about 2 – 5 million/mL. All samples must be filtered (minimum 40 micron pore) before sorting.
2. **Bring extra pre-sort buffer.** Buffer in which to bring your cells, although some variation is acceptable:
 - a. *Optimal* - BD Biosciences Pre-sort buffer (cat#563503)
 - b. 1 x Phosphate Buffered Saline (PBS) +1% Fetal Calf Serum (Heat inactivated) or 1% Albumin (for sticky cells)
 - c. Hanks Balanced Salt Solution (HBSS) (Ca²/ Mg² Free) 25 mM HEPES pH7.0
1% Fetal Calf Serum (Heat inactivated) or 1% Albumin (for sticky cells)
3. For samples with a high dead cell concentration or primary tissue preps: Dead cells can release their DNA into sorting media which in turn can cause cells to clump together.

Strategies to reduce cell clumping

- a. BD Biosciences pre-sort buffer
 - b. Include 50% Accumax (Invitrogen) in sort buffer
 - c. Treat cells for 15-30 minutes in a sterile solution of 100 ug/mL DNase I and 5 mM MgCl₂ in HBSS at room temp. Wash the cells 1x in HBSS containing 5mM MgCl₂. Resuspend the cells in HBSS containing 25-50 ug/mL DNase I, plus at least 1mM MgCl₂ prior to and during the sort. (5mM MgCl₂ is optimal)
4. For bulk collection, **bring media with a high FBS content (20-50%)** to help cells recover from the sorting process.
 - a. For direct plate collection, bring plates loaded with the culture media appropriate for the sorted cells.
 - b. For post sort mRNA analysis, cells can be sorted directly into RNA isolation reagents, such as Trizol-LS, RTL buffer, RNA Later, etc. Please ask for more details.